

General Protein Swabs

Emport
LLC

Swabbing surfaces for general protein is a quick, cost-effective way to monitor allergen cleanliness. Protein swabs are often used to verify day-to-day cleaning activities, providing presence/absence results with little-to-no equipment or waiting time. The swabs will react to any form of protein, whether allergenic or not, with a color change. They are generally less sensitive and specific than allergen-specific rapid tests, and so they are not suitable for program validation.

Hygiene

AllerSnap

- 100 swabs/case
- High sensitivity
- 5-15 minutes to results
- Long-handled snap bulb design
- Requires incubator

PRO-Clean

- 100 swabs/case
- 1-10 minutes to results
- Long-handled snap bulb design
- No equipment required

Gold Standard Diagnostics

Path-Chek Protein

- 50 double-sided swabs/case
- Instant results
- Short swab design
- No equipment required

Bio-Check UK

Protein Swab

- 25 swabs/case
- Instant results
- Long-handled snap bulb design
- No equipment required
- Refrigerated storage

FOR MORE INFORMATION
EmportLLC.com | info@emportllc.com
412.447.1888 | 866.509.4482 toll free



For Details

