

: OnSite® Gluten

The OnSite® Gluten Test kit is designed for any facility that needs rapid, accurate gluten detection: manufacturers, commercial kitchens, food analysis laboratories, and others. The kit contains everything needed to detect traces of gluten in ingredients, finished products, and the working environment.

: OnSite® Gluten

- **Adjustable LOD:** set product sensitivity to 5, 10, or 20 mg/kg (ppm)
- **Use** for foods, beverages, rinse water and surfaces
- **No** lab equipment required; swabs included in the box
- **Results** in <20 minutes
- **2D4 antibody** technology detects gluten from both raw/native and processed (deamidated) foods
- **Matrices validated** according to AOAC PTM protocols:
Stainless steel surface, spice mix, rice flour, oat flour, baked bread.

Additional testing kits include OnSite® III Allergen and OnSite® Allergen Multiplex Tree Nut.

FOR MORE INFORMATION
EmporLLC.com | info@emporllc.com
412.447.1888 | 866.509.4482 toll free



Empor
LLC