


AlerTox® Sticks

**Test ingredients, finished foods & surfaces for
Nuts, Crustacean, Egg, Fish, Soy & More!**

AlerTox Sticks kits are specially designed to help food manufacturers monitor levels of common allergenic proteins at every step of production and in a variety of matrices. Kits are user-friendly, accurate, and robust: a must-have for any manufacturer or commercial kitchen that wants to ensure food allergen safety.

Allergen	Strips/kit	LOD*
AlerTox Sticks Almond	10 strips	20 ppm
AlerTox Sticks Crustacean	10 strips	10 ppm
AlerTox Sticks Egg	10 strips	1.25 ppm
AlerTox Sticks Fish	10 strips	5 ppm
AlerTox Sticks Hazelnut	10 strips	20 ppm
AlerTox Sticks Mustard	10 strips	2 ppm
AlerTox Sticks Peanut	10 strips	Varies: ask for details
AlerTox Sticks Soy	10 strips	10 ppm
AlerTox Sticks Total Milk	10 strips	2.5ppm



AlerTox Sticks are:

- ◆ **Quick:** A reliable answer in less than 10 minutes
- ◆ **Easy:** No special equipment required; everything is included in the kit
- ◆ **Reliable:** Shelf-stable and long-lasting
- ◆ **Cost efficient:** competitive pricing and volume discounts available
- ◆ **ELISA kits** also available for 20+ allergens

*LOD may vary in some matrices; see manual for details. AlerTox Sticks are designed to detect trace contamination. At very high concentrations of allergen, a hook effect may be present. Inquire for details.

Looking for more test kits? Also available:



GlutenTox

→ Rapid onsite kits for gluten detection in ingredients and surfaces

**5-in-1
Water Test**

→ Quickly test water for pH, nitrites, nitrates, & more

oleotest
TRUST IT

→ Test for polar compounds to keep your oil fresh

Ready to order? Contact Emport LLC:

orders@emportllc.com · 412-447-1888 · 866-509-4482 · EmportLLC.com