

General Instructions for Use Qualitative Gluten PT Standard Set

DO NOT OPEN VIAL UNTIL INSTRUCTIONS HAVE BEEN REVIEWED.

Sample Description

- This PT sample set is designed for the qualitative determination of gluten allergen.
- The PT sample is provided as a 3 sample set. Each vial contains approximately 10 grams of material.
- At least one vial on the set will be positive for the allergen of interest.
- The PT material has been evaluated for compatibility with multiple Elisa type assays.
- Store the PT sample at 15-30°C until you are ready to analyze.

Preparation Instructions

- Retrieve samples and allow them to equilibrate to room temperature.
- Mix each sample by gently shaking or rolling the vial.
- Subsample appropriate aliquot from each vial for your analytical method.
- **Report “1” for present and “0” for absent for the allergen of interest for each sample.**