



February 2, 2024 — Emport LLC offers new OnSite® Allergen Multiplex Tree Nut test kit for food manufacturing facilities and commercial kitchens.

The OnSite® Allergen Multiplex Tree Nut test kit offers the food and hospitality industries an opportunity to quickly and easily monitor for multiple common tree nuts in rinsewater and on surfaces.

Emport LLC is pleased to stock the first kit from OnSite® Food Safety, titled OnSite® Allergen Multiplex Tree Nut. This rapid lateral flow test allows users to simultaneously detect almond, cashew, hazelnut, pecan, pistachio, and walnut contamination on surfaces and in rinsewater. It is a natural fit within Emport LLC's catalog of easy, accurate allergen management test kits.

Careful planning and daily vigilance are necessary to safeguard the millions of people worldwide who live with food allergies. In the United States and abroad, government regulations frequently require careful management of common allergens including tree nuts. In addition, voluntary but popular GFSI food safety schemes require comprehensive allergen management. The OnSite® Allergen Multiplex Tree Nut kit is an important daily tool for food companies and commercial kitchens alike.

The OnSite® Allergen Multiplex Tree Nut rapid test requires no equipment and provides results in less than 20 minutes. The test is sensitive enough to detect as little as 1-5 mg/kg (ppm) of tree nut protein in rinsewater or 10µg/100 cm² on surfaces. If any or all of the targeted tree nut proteins are present in sufficient quantity, the test cassette will display a red test line alongside the internal control line.

OnSite® Allergen Multiplex Tree Nut is an affordable part of any allergen management plan. It can be paired with a variety of single-nut rapid tests and/or general protein swabs from Emport LLC's catalog, as well as with third party laboratory testing.

For more information, including purchasing information, visit emportllc.com.

About OnSite®

OnSite® is a new line of food safety test kits produced by Microbiologique. OnSite® tests are simple, rapid methods designed for laboratory and industry use.

About Emport LLC

Since 2011, Emport LLC has operated with a simple mission: More Safe Food, More Happy People. The company distributes a growing assortment of food safety supplies, with a focus on proteins including gluten, allergens, and meat. Emport LLC's team members specialize in finding the right products for each customer, each time.