



## Path-Chek® Hygiene Protein

Path-Chek® Hygiene Protein is an inexpensive and rapid chemical based test for the detection of invisible food residues on food manufacturing surfaces after routine cleaning and sanitising procedures have been performed. The test provides an instant result without the need for laboratory facilities or incubation.

Path-Chek® Hygiene Protein does not detect bacteria, so tests for specific food borne pathogens should always be performed, but protein and food residues detected by the test indicate that bacteria may also be present on the surface.

Path-Chek® Hygiene Protein swabs are economical to use as one swab has two swab heads and therefore two tests per swab (opened swabs can be stored in a zip lock bag for up to 7 days at 2-8°C).

### Key Benefits



**Results in 5 seconds**



**Simple to use**



**Cost effective -  
2 tests per swab**



**No laboratory facilities  
required**